## New methods for iron isotope measurements in ancient sedimentary rocks Aleksandra Galic<sup>1</sup>, Paul R.D. Mason<sup>1</sup> and Pieter Z. Vroon<sup>2</sup> <sup>1</sup>Department of Earth Sciences, Utrecht University, Budapestlaan 4, 3584 CD Utrecht, <sup>2</sup>Faculty of Earth and Life Sciences,

Vrije Universiteit, De Boelelaan 1085, 1081 HV Amsterdam

Fe isotopes are likely to be an important tracer for finding traces of possible life during future missions to other planets in our solar system

- Current Mullticollector ICP-MS-based measurements techniques are prone to substantial bias due to inadequate sample preparation.
- Fe isotopes vary over a limited range in nature (order of 3.5‰), and highly precise and accurate isotopic measurements are thus essential

## **RESULTS:**

Four different methods were performed to test 100% recovery of iron:

1. Elution with 6.0M HCl:

#### ~2% Fe lost in prefraction

- 2. Elution as in 1) with 6.0M HCl +0.001%  $H_2O_2$
- 3. Elution as above but with different resin particle size
- 4. Elution of the prefraction was performed testing different chemical
- Here we aim to obtain a 100% Fe yield after column extractions and completely remove all possible interfering elements to give accurate results

## **MATERIALS & METHODS:**

- We explored a range of previously published techniques (Dauphas et al. 2009, Dauphas and Rouxel, 2006 and Schoenberg and von Blanckenburg, 2005) for Fe extraction from Fe-rich rocks including banded iron formations, shales and pyrite
- ♦ Two different resin types were tested: AG1-X8 200-400 and AG1-X8 100-200 mesh Cl<sup>-</sup> form anion exchange resin.

Artificial solution

Rock reference material BHVO-2 (basalt)

Mineral sample

### (6.0M HCl + 0.01M HClO<sub>4</sub>) Elution with perclorate; Ti follows Fe:



**Reference material test:** 

Twodifferentelutiontechniquesweretestedandtwocolumnextractionsperformed:



	→ 8.63% Fe	Pyrito
10 ppm Fe	MAG-1 (fine grained clayey mud)	i ynte
2ppm Ca	→ 4.75 wt% Fe	Particle dimension: 60-120
2ppm Cr	ECSC 681-1 (iron ore)	μm
2ppm Ni	→ 33.21wt % Fe	~ 45 wt% Fe
2ppm Ti		

♦ Two different oxidizing agents ( $H_20_2$  and  $HClO_4$ ) were used to keep Fe in its Fe<sup>3+</sup> form in order to obtain a 100% recovery of iron.





# HCl + 0.001%H<sub>2</sub>O<sub>2</sub> HCl+ 0.01M HClO<sub>4</sub>

Fe yield ~100% Ti is completely eliminated after the 2<sup>nd</sup> column extraction

**Pyrite test:** 

Pyrite samples are dissolved in:

•10 drops HF conc., 6 drops HNO3 conc. and 1 drop HClO4 conc.
•Prefraction eluted with 6.0M HCl+0.01M HClO<sub>4</sub>
•Fe fraction eluted with 0.4 M HCl



Fe yield ~ 100%, no S

*Iron isotpes will be measured on materials obained from a International Continental Drilling Program (ICDP) in Barbeton Greenstone Belt in South Africa* 

#### References cited: **Dauphas et al.**, Chemical Geology 267 (2009), **Dauphas and Rouxel**, Mass Spectrometry Reviews 25 (2006), **Schoenberg and von Blanckenburg**, International Journal of Mass Spectrometry 25 (2005)



#### ✤ ~ 100%Fe yield in all tests

- ✤ HCl only → risk of losing some Fe in prefraction
- HClO<sub>4</sub> and H<sub>2</sub>O<sub>2</sub> agents result both in 100% yields. However HClO<sub>4</sub> is preferred because it is more stable during the evaporation of solutions in subsequent steps of the method.
- **\*** Ti follows Fe chemistry but gets eliminated in 2<sup>nd</sup> extraction
- The pyrite test resulted in good yields without sulpher in the final iron solution
- The collection of the Fe from the columns can be done without HClO<sub>4</sub> for pyrite samples